



August 1, 2009

To Whom It May Concern:

Testing performed by Micrus Endovascular Corporation demonstrates that Micrus® Microcoils are MRI compatible and show no appreciable temperature rise, movement, or artifact when exposed to Magnetic Resonance Imaging (MRI). The Micrus Microcoil Systems include CERECYTE® and Stretch Resistant Platinum Microcoils, as listed below:

- MICRUSPHERE® Microcoils
- PRESIDIO® Microcoils
- CASHMERE™ Microcoils
- DELTAPAQ™ Microcoils
- DELTAPLUSH™ Microcoils
- ULTIPAQ® Microcoils
- HELIPAQ® Microcoils
- INTERPAQ® Microcoils.

Diagnostic MRI compatibility with Micrus Microcoils was tested at field strengths up to and including 3.0 Tesla.

If you require additional information, please do not hesitate to refer to the instructions for use, which is included with each microcoil system, or contact our Customer Service department at 1-888-550-4120 ext. 443 in the US, or 1-408-433-1443 outside the US.

Sincerely,

A handwritten signature in black ink, appearing to read 'R. Michael Crompton', written over a large, stylized flourish that extends to the right.

R. Michael Crompton
Vice President, Regulatory/
Clinical Affairs & Quality

P/N 50597-01 Rev E